

Bead View Report Generator

Refer to the users manual pages listed for more in depth information.

Open the Bead View Software (Page 5)

1. Click "File"
2. Select "Open"
 - a. Desktop-My Computer-Local Disk-My Batches
3. Select "Output"
4. Click "Open"

Plate Icon (Pages 6-7)

1. In the Drop Down Menu-Select Type
2. Select well
3. Enter the highest concentration as S7
4. Continue down to S1
5. Enter SB in the blank

Samples Icon (Page 9)

1. Fill in the chart
2. Mark your duplicates

Kit Set-Up (Page 10)

1. Analyte Tab-select Background Subtraction-yes
2. Standards-
 - a. Under S7 enter the full concentration of your standard
 - b. Continue through S1- the most dilute
3. Recommended Formula-Select Best Fit

Analysis Icon (Page 12)

1. Click "+*" calculates all analytes

Replicates (Page 13)

1. Click "Generate Report"